


<Luciferase assay>

0. Promoterに対する転写因子、転写調節因子の活性をLuciferaseをReporterとして測定する。



1. 準備

Lysis buffer (1x)

2. Procedure

Day 1

細胞をまく { $1 \sim 4 \times 10^5$ cells / well (12 or 24 well plate), duplicateで行う}

293T : 2×10^5 cells / well / 12 well plate or 1×10^5 cells / well / 24 well plate

HepG2 : 4×10^5 cells / well / 12 well plate

Mv1Lu : 2×10^5 cells / well / 12 well plate

R4-2 : 2×10^5 cells / well / 12 well plate or 1×10^5 cells / well / 24 well plate

Day 2

Transfection

Reporter (-luc)

pRL-TK, pRL- CMV, β -gal (CH110) (補正用)

Day 3

Cells (12 or 24 well plate)

↓ PBS wash x 1 (293Tの場合は省略)

↓ add 150 μ l of Lysis buffer

↓ rotate 20 min, RT (この時発光試薬、2 x β -gal solu.を溶かしておく)

↓ to microtube

↓ 15000 rpm, 10 min, 4°C

sup

<Luciferase assay>

15 μ l

↓ to 96 well plate (白色)

発光試薬を機械に設置し測定①

< β -gal assay (Transfection効率の補正) >

30 μ l

↓ to 96 well plate (透明)

↓ add 30 μ l of 2 x β -gal solu.

↓ 37°C, 0.5h ~ O/N (液が黄色になるまで)

2Fのplate readerで測定②

計算

① / ②を計算 (s.dをつける)