

Preparation of plasmid DNA (Mini-Prep)

(用途)

インサートのチェックの時に用いる。

(とったPlasmidをシーケンスするときはPCR purification kitで精製)

(Day1; **Bacteria culture**)

Colony (LB / Amp plate) → 2 ml of LB (or 2 x YT) / Amp (50-100 $\mu\text{g/ml}$)
in 15 ml tube

↓ Culture O/N (~ 18 h)

(Day2; **Mini-prep**)

Bacteria in 15 ml tubes

↓ to 1.5 ml tube
↓ 15,000 rpm, 1 min, 4°C
ppt
↓ add 200 μl of Sol. I
↓ suspend by vortex
↓ add 200 μl of Sol. II (Lysis buffer)
↓ invert 7 times (clear cell lysate)
↓ add 200 μl of Sol. III (Neutralization buffer)
↓ invert 5 times (become clouded)
↓ Add 10 μl of chloroform
↓ 15,000 rpm, 10 min, 4°C
sup (~ 550 μl)
↓ add 550 μl of isopropanol
↓ R.T 5 min
↓ 15,000 rpm, 10 min, 4°C
ppt
↓ add 500 μl of 70% EtOH
↓ vortex
↓ 15,000 rpm, 5 min, 4°C
ppt
↓ air dry
↓ dissolve in 40 μl of TE

Sol. I
50 mM Tris-HCl (pH 7.5)
10 mM EDTA
50-100 $\mu\text{g/ml}$ RNase A

Sol. II
0.2M NaOH
1% SDS

Sol. III
1.32M KOAc (pH 4.8)

Quantification (1: 20) & cut check (5~10 μl)