

## <Sphere formation assay>

### 1. Preparation

#### i) Sphere formation medium

<for glioma cells (standard)>

DMEM	25 ml (with antibiotics)
F12 medium	25 ml (with antibiotics)
20 ng/ml of EGF	10 $\mu$ l of stock (100 $\mu$ g/ml)
20 ng/ml of bFGF	10 $\mu$ l of stock (100 $\mu$ g/ml)
B27 supplement	1 ml

Cells: U251,  
NMuMG

<for epithelial tumor cells>

DMEM	25 ml (with antibiotics)
F12 medium	25 ml (with antibiotics)
20 ng/ml of EGF	10 $\mu$ l of stock (100 $\mu$ g/ml)
20 ng/ml of bFGF	10 $\mu$ l of stock (100 $\mu$ g/ml)
B27 supplement	0.5 ml

Cells: TE13, HaCaT-Ras<sup>G12V</sup>  
154, 234

#### ii) Dishes

Epithelial tumor cells: Ultra-Low Attachment 10 cm dish (Corning 3262)  
Ultra-Low Attachment 6 well plate (Corning 3471)

Glioma cells: Non coated dish (for bacterial culture)

### 2. Procedure

Cells

- ↓ Suspend in sphere formation medium and count
- ↓ Seed in 10 cm ( $5 \times 10^4$  cells) or  
6 well plate ( $5 \times 10^3$  cells)
- ↓ 5 ~ 14 days (depend on cell type)
- ↓ Fix with formaldehyde (final 4%)

Microscope OR count visible spheres

If spheres are fused with each other, please consider the addition of methylcellulose (final 1%) to medium.